

LT-IWDO 1"/0.68 low temperature compatible aspheric objective for attoCFM

I,
1021992

Technical Specifications

| | |
|--------------------------------------|--|
| Optical Data | |
| clear aperture [mm] | 5 |
| focal length [mm] | 3.1 |
| numerical aperture(NA) | 0.68 |
| working distance | 1.5 mm |
| Spectral Performance | |
| AR coating (> 80% transmission) [nm] | A : 350 .. 700 B : 600 .. 1050 C : 1050 .. 1700 |
| Compatibility | |
| environment | high magnetic fields high vacuum low temperature |
| compatible setups | attoCFM I, attoCFM IV, attoDRY800 |
| Size and Dimensions | |
| length | 51 |